

DIAGNOSTIC WAX UP CHECK LIST

To help us serve you better KYORI DENTAL LAB has put together a check list of the items needed to complete a diagnostic wax up. These diagnostic tools should be obtained during the diagnostic appointment a couple of weeks before the prepping appointment.

1. Photos of pre-op teeth including: Nose to chin natural smile, nose to chin retracted, 1:2 retracted, profile of natural smile. (Please use translucent cheek retractors.)
2. Bite registration including CR/CO - (Centric Relation/Centric Occlusion)
3. Stick-bite - vertical and horizontal planes (canter and midline)
4. Photo of stick bite in place
5. Impression of upper and lower pre-op. (Palate and Gum tissue is as important as the capturing of a margin in the restoration.)
6. Indicate the desired length of final restorations (with temp Impression)
7. Vertical measurements from CEJ to CEJ - (Cemento-Enamel Junction)
8. Detailed lab prescription, expressing goals of patient, thoroughly filed out.

HERE ARE A COUPLE OF HINTS WHEN SENDING A CASE:

1. Send models and polyvinyl Impressions.
2. Label all photos with the doctor's name, patient's name and date.
3. Make sure the models are securely packed or they will brake.
4. Make sure the Bite Registration and Stickbite are securely packed. If the Stickbite moves it becomes useless.

This checklist is necessary to help us serve you and your patients better.

KYORI DENTAL LAB thanks you.